Reactions involved:

Page No. ____18_

Experiment 4

Aim:

To determine the amount of fermanent, temporary and total hardness of given water sample using standard N/50 EDTA solution.

Affaratus required:

Measuring cylinder, fifette, wash bottle, standard flask, burette, beaker, conical flask.

Reagents required:

Water, EDTA, Eriochrome black T, buffer

Theory: Water that doesn't produce lather with soap solutions nother produces white precipitate or scum is known as hard water. Hardness in water is due to the presence of dissolved salts of calcium and magnesium in the form of bicarbonate, chloride & sulphate in water It is unfit for drinking, bathing, washing and it also forms scales in boilers. Hence, it is necessary to estimate the amount of hardness producing substance present in water sample Once it is estimated, the amount of chemicals required for the treatment of water can be calculated. There are two types of hardness fresent in water: temporary hardness & fermanent hardness. Temporary hardness can be removed by boiling the water or by the addition of lime (or calcium hydroxide). The estimation of hardness is based on complexometric titration



Teacher's Signature: ____

Observations: Titration of given water sample against EDTA solution for total hardness

S.no.	Initial burette readings = V.	Final burette readings = Ve	Volume of given NazEDTA sol used = (Vz-V1)
1.	0.0 ml	3.0 ml	3.0 ml
2.	3.0 ml	6.0 ml	3.0 ml
3.	6.0 ml	9.0 ml	3.0 ml

Titration of given water sample against EDTA solution for fermament hardness

s.no.	Initial burette readings = Vi	Final burette readings = V2	Vol ^m of standardized EDTA sol" wed = (V ₂ -V ₁)
1.	0.0 ml	2.5 ml	2.5 ml
2.	2.5 ml	50ml	2.5 ml
3.	5.0 ml	7.5 ml	2.5ml

Calculations:

Total hardness:

$$N_{EDTA} = \frac{N}{50}$$

Number =
$$\frac{N_{\text{EDTA}} \times V_{\text{EDTA}}}{V_{\text{uniter}}} = \frac{1}{50} \times \frac{3}{20} = \frac{3}{1000} \text{ N}$$

	Date
Expt. No./ Name:	Page No J9
Hardness of water is determine of ethylene diamine tetra ace agent. Since EDTA is insoluble is taken for this experiment. I bonds with a metal ion.	ed by titrating with a standard solution lic acid (EDTA) which is a complexing. in water, the disodium salt of EDTA DTA can form four or six coordination
Procedure:	Paramet Langue ?
Total hardness of the water sa	mble:
in Pibetted out 20 ml of water sa	mble: mble into a 250 ml conical flask Added uffer & I drop of EBT indicator You got a
2-3 ml of SH 10 ammonical by	uller & I drop of EBT indicator You got a
(i) Titrated water sample with E	DTA solution until the colour changed from
wine nod to how Recorded the	roume
Refeated this 2 more times till	you got concordant readings.
The state of the part of the p	and of together was wear Improved
2. Permanent hardness:	10. + 05. 01 A. C. +0 L.: (1.)
is Transferred a 100 ml water &	ample into 250 ml beater trendy soula
the water for 30 minutes. Allo	ample into 250 ml beaker Grently boiled water to cool to room
T- 1410-4	
Filtered the water directly in	to a clean 100 ml volumetric flask. Filled ark with distilled water.
the volumetric flash to the mi	sample into a flack & added 2-3 ml
Pifetted 20 ml of the julioned	1 dual of ERT indicator.
of BH 10 ammonical buffer &	th EDTA rolution until the colour changed
	A J TIM A LYBY LL MOU
w Repeated this 2 more times til	I you got concordant readings.
w Repeated this 2 more willes we	3
	Teacher's Signature :
	leacher 3 digitature

Total hardness with respect to $CaCO_3 = N_{water} \times Eq. wt. of <math>CaCO_3 = \frac{3}{1000} \times 50 = 0.15 \text{ g L}^{-1}$

Total hardness in frm = Nwater × Eq. wt. of CaCO3 × 1000 = 0.15 × 1000 = 150 frm

Permanent hardness:

Varater = 20 ml ,, VEDTA = 2.5 ml

 $N_{\text{EDTA}} = \frac{N}{50}$

Nwater × Vunter = $N_{EDTA} \times V_{EDTA}$ Number = $\frac{N_{EDTA} \times V_{EDTA}}{V_{water}} = \frac{1}{50} \times \frac{2.5}{20} = \frac{2.5}{1000}$

Permanent hardness with respect to $CaCO_3 = N$ water × Eq. wt. of $CaCO_3 = \frac{2.5}{1000} \times 50 = 0.125 \text{ g L}^{-1}$

Permanent hardness in ppm = Nwater × Eq. wt. of Ca CO3 × 1000 = 0.125 × 1000 = 125 ppm

Temporary hardness: Temporary hardness in from = Total hardness - Permanent hardness = 150 - 125 = 25 from

	Date
Expt. No./ Name:	Page No20
Result:	
in Total hardness = 150 ppm	
in Permanent hardness = 125 pm	
Temporary hardness = Total har	dress - Permanent hardness = 25 fpm
Precautions:	
is Before starting the experiment,	the glass affaratus must be ferfectly
in Always read lower meniscus of	solution level in purette
Near end point, add EDTA volu	rolution level in burette. tion drop wise and after addition of each
drop, see the colour change	ion from fifette. Just tap the tip of the
w Do not blow last drop of rolut	ion from hibette. Just tap the tip of the
sifethe to the walls of the flash	0 11
	Teacher's Signature :